

Alcohol Testing Program - Instrument Processing Sheet

Agency: Jacksonville SO

Instrument Serial Number: 80-001743

Date In: 8/27/2026

DI Completion Date: 09/01/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: IAS Date: 8/27/2026		Quality Checks By: IAS Date: 8/27/2026		Flow Adjustment By: _____																
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☒ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 220 ☒ Flow Verification (L/s) Flow Column #: ATP102 32 mm .156 (.139-.169) 36 mm .171 (.156-.190) 53 mm .242 (.228-.278) 103 mm .519 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 28427 Gauge: 1010 Instrument: 1010 ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm (.139-.169) 36 mm (.156-.190) 53 mm (.228-.278) 103 mm (.447-.547)																
		<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>202603A 3/3/2028</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>202603B 3/3/2028</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>202603D 3/4/2028</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>23725080A1 10/5/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202603A 3/3/2028	0.080	MP5089	202603B 3/3/2028	0.200	MP5090	202603D 3/4/2028	0.080 DGS	N/A	23725080A1 10/5/2027	Maintenance By: _____ Date: _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other:	
Simulator	Serial #	Lot#/Exp																		
0.050	MP5088	202603A 3/3/2028																		
0.080	MP5089	202603B 3/3/2028																		
0.200	MP5090	202603D 3/4/2028																		
0.080 DGS	N/A	23725080A1 10/5/2027																		
Optical Bench Adjustment By: _____		Department Inspection By: SLH																		
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: 28427																		
Simulator	Serial #	Lot #	Expiration	Gauge: 1012 Instrument: 1012 Mouth Alcohol Solution Lot #: 2025-D Exp: 09/25/2027 Acetone Stock Solution Lot #: 2025-B Exp: 09/22/2027																
0.000		N/A	N/A																	
0.040																				
0.100																				
0.200																				
0.300																				
0.080 DGS	N/A																			
☐ Post Optical Bench Adjustment Stability Checks																				
Simulator	Serial #	Lot #	Expiration																	
0.050																				
0.080																				
0.200																				
0.080 DGS	N/A																			
Gauge ID #: _____																				
Gauge: _____ Instrument: _____																				
Notes/Suggested Service: Instrument was able to read the DGS psi. SLH 9/1/2026				☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment ☒ Average Calculations																
				☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☐ Other:																
				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																
				Digitally signed by Wen-Chi Pierson Date: 2026.09.09 08:47:50 -04'00' Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.09.10 07:23:53 -04'00'																
				Tech Review Admin Review																

Stability Checks

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.205 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L 50.003 g/210L of Wet
POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Test g/210L Time Air Blank 0.000 15:03 Control Test 0.050 15:04 Air Blank 0.000 15:05 Control Test 0.050 15:06 Air Blank 0.000 15:07 Control Test 0.050 15:07 Control Test Stats Average 0.0500 Std Dev 0.0000 Rel. Std Dev(%) 0.0000	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Test g/210L Time Air Blank 0.000 15:40 Control Test 0.079 15:41 Air Blank 0.000 15:41 Control Test 0.078 15:42 Air Blank 0.000 15:42 Control Test 0.079 15:43 Air Blank 0.000 15:43 Control Test Stats Average 0.0787 Std Dev 0.0006 Rel. Std Dev(%) 0.7335	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Test g/210L Time Air Blank 0.000 15:52 Control Test 0.199 15:52 Air Blank 0.000 15:53 Control Test 0.198 15:54 Air Blank 0.000 15:54 Control Test 0.198 15:55 Air Blank 0.000 15:55 Control Test Stats Average 0.1983 Std Dev 0.0006 Rel. Std Dev(%) 0.2911	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Test g/210L Time Air Blank 0.000 15:57 Control Test 0.081 15:58 Air Blank 0.000 15:58 Control Test 0.080 15:59 Air Blank 0.000 15:59 Control Test 0.081 15:59 Air Blank 0.000 16:00 Control Test Stats Average 0.0807 Std Dev 0.0006 Rel. Std Dev(%) 0.7157
POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Operator's Signature	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Operator's Signature	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Operator's Signature	POCONO, PA 18342 Intelligence - Ricordi Analyser Model 8100 08/27/2026 Software: 6100.27 Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: JACKSONVILLE SO
Time of Inspection: 17:56

Date of Inspection: 09/01/2026

Serial Number: 80-001743
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:28424080A3 Exp: 11/05/2026
0.000	0.050	0.078	0.197	0.081
0.000	0.049	0.078	0.197	0.081
0.000	0.050	0.078	0.197	0.081
0.000	0.049	0.078	0.197	0.081
0.000	0.050	0.079	0.197	0.081
0.000	0.049	0.079	0.197	0.081
0.000	0.050	0.079	0.197	0.080
0.000	0.050	0.079	0.197	0.081
0.000	0.050	0.079	0.197	0.080
0.000	0.049	0.079	0.198	0.081


Standard Deviations	0.0005	0.0005	0.0003	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


LEANDRA HIGGINBOTHAM
Signature and Printed Name

09/01/2026
Date

Florida Department of Law Enforcement
Alcohol Testing Program - Average Calculations for Calibration Certificate

Calculations By: SLH

Serial Number: 80-001743

Date: 9/1/2026

Solution:

0.05 g/210 L Test: 0.496 / 10 = 0.050 g/210L
(SUM) (AVERAGE)

0.08 g/210L Test: 0.786 / 10 = 0.079 g/210L
(SUM) (AVERAGE)

0.20 g/210L Test: 1.971 / 10 = 0.197 g/210L
(SUM) (AVERAGE)

0.08 g/210L Dry Gas Std Test: 0.808 / 10 = 0.081 g/210L
(SUM) (AVERAGE)

Comments:

Average = (Sum of observations) / (Total number of observations)

All results are reported in g/210L



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number **80-001743**, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	80-001743	Average of y ± Uncertainty*
Owning Agency:	JACKSONVILLE SO	0.050 g/ 210 L 0.004
Calibration Date:	09/01/2026	0.079 g/ 210 L 0.004
Calibration Time:	17:56	0.197 g/ 210 L 0.008
		0.081 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

S. LeAndra
Higginbotham

Digitally signed by S. LeAndra
Higginbotham
Date: 2026.09.04 13:44:10 -04'00'

09/01/2026

Date

LEANDRA HIGGINBOTHAM,

Department Inspector

FDLE/ATP Form 69 August 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality